



Goat anti-Mouse IgG (H+L), FITC AffiniPure

Catalog	PMK-014-093M	PMK-014-093S	Tel :027-87888039
Quantity	100μL	1ml	E-mail:postmaster@biopmk.com
			Web: www.biopmk.com

For research use only.

Applications	Species Cross-Reactivity	Molecular Weight	Isotype
Flow Cyt, ICC/IF, IHC	Mouse	N/A	IgG1

Storage Buffer & Condition: 0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6

Store at -20°C. Do not aliquot the antibody.

Stabilizer: 15mg/ml Bovine Serum Albumin (IgG-Free, Protease-Free)

Preservative: 0.05% Sodium Azide

Suggested Dilution Range: 1:100-1:800 for most applications

Optimal dilutions should be determined by the end user.

Fluorophore: Fluorescein-5-isothiocyanate (FITC)
Amax=492nm; Emax=520nm

Purity: The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

Antibody Specificity: Based on immunoelectrophoresis and/or ELISA, the antibody reacts with whole molecule mouse IgG. It also reacts with the light chains of other mouse immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins.

The antibody may cross-react with immunoglobulins from other species.

Expiration date: Store at -20°C. Avoid repeated freeze / thaw cycles. The expiration date may be extended if test results are acceptable for the intended use.

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.